


# SmartSoc BYOD

## What is BYOD?


- Bring Your Own Device (BYOD) is a SmartSoc update option now being offered to customers to turn their iPhone 11, iPhone 11 Pro, or iPhone 11 Pro Max into a SmartSoc scanning device.
- Currently, only iPhone 11, iPhone 11 Pro, and iPhone Pro Max models are supported for this BYOD option.  
- Ref. **SmartSoc BYOD - Approved Device Models (SL00266 Rev A)**
- If a BYOD approved device cannot be obtained, your SmartSoc system can still be upgraded by purchasing a device through Orthomerica.  
- Ref. **SmartSoc by Orthomerica (SL00264 Rev A)**

## Upgrade a device through BYOD: How To


1. Acquire an approved BYOD mobile device. - Ref. **SmartSoc BYOD - Approved Device Models (SL00266 Rev A)**
  - a. **If the device is brand new**, it will need to be setup for personal use before continuing.  
- Ref. **SmartSoc BYOD – New Device Setup (SL00267 Rev A)**
  - b. **If the device is currently operated for personal use**, continue to download - **valid Apple ID is required.**
2. Download the “**CurveCapture O&P**” (v4.06) scanning application:
  - a. Launch the iOS “**App Store**” on your iPhone




b. Select the "Search" tab




c. Search "CurveCapture O&P" in the search bar



d. Select the mobile app **“CurveCapture O&P”** (v4.06)



e. Download **“CurveCapture O&P”** (v4.06)



- 3. Complete the following provided forms:
  - a. SmartSoc Amendment 3
    - Ref. **SmartSoc Amendment (FM013.09 Rev 3)**
  - b. SmartSoc Order Form for BYOD
    - Ref. **SmartSoc Order Form (BYOD) (FM002.101 Rev A)**

**SmartSoc Order Form (BYOD) – How To**

**SmartSoc BYOD Order Form**

**REQUIRED**  
Company

**REQUIRED**  
CPI Account Number

**REQUIRED**  
Billing Address

**REQUIRED**   **REQUIRED**   **REQUIRED**  
City                      State                      Zip

**REQUIRED**  
CPI Sales Rep

**REQUIRED**  
PO#




Practitioners must pass Orthomerica's STARband® Cranial Course (SCC) certification course or equivalent in order to purchase Cranial Remolding Orthoses. Please visit [www.orthomerica.com/education](http://www.orthomerica.com/education) for more information.

| PRACTITIONER   |                          | DEVICE                                    |
|--|--------------------------|---|
| <b>REQUIRED</b><br>Name  | <b>REQUIRED</b><br>Phone | <b>REQUIRED - Ref. 1</b><br>Apple ID      |
| <b>REQUIRED</b><br>Email   |                          | <b>REQUIRED - Ref. 2</b><br>Model Name    |
| <input type="checkbox"/> Cranial Module <input type="checkbox"/> Prosthetic Module <b>Required for Cranial Module only</b> |                          | <b>REQUIRED - Ref. 3</b><br>Model Number  |
| <b>REQUIRED</b><br>Certification Number  |                          | <b>REQUIRED - Ref. 4</b><br>Serial Number |
| <b>REQUIRED</b><br>Shipping Address  |                          | <b>REQUIRED - Ref. 5</b><br>Device Serial |
| <b>REQUIRED</b><br>City  | <b>REQUIRED</b><br>State |   |
| <b>REQUIRED</b><br>Zip   |                          |   |
| <b>REQUIRED</b><br>C-Fab Company   |                          |   |
| <b>REQUIRED</b><br>C-Fab Contact   |                          |   |
| <b>REQUIRED</b><br>C-Fab Email (scans will be submitted to this address)   |                          |   |


▼ **If a scan must be rushed, Orthomerica must be notified via email or phone call to expedite the order** ▼

## 1 How to find “Apple ID”:

1. Open “Settings”
2. Select “Apple ID” (first option)




3. Record Apple ID email address:




## 2, 3, 4 How to find “Model Name”, “Model Number” and “Serial Number”:

1. Open “Settings”
2. Select “General”
3. Select “About”
4. Record information:



## 5 How to find “Device Serial”:

1. Open “CurveCapture O&P”
2. Select :



**Note:** For additional assistance, please view the video guide.


– Ref. **Video - How to Setup Your BYOD Device (V08.002.05 Rev A)**

4. Submit the completed forms via email to the following: [custserv@orthomerica.com](mailto:custserv@orthomerica.com) .
5. Upon receiving, Orthomerica will process the registration and advise to proceed with activating your system.
  - a. User will receive the activation code to activate the device through the *CurveCapture O&P* mobile app.
  - b. Upon activation, the *eBrace Global* web-portal and the *CurveCapture O&P* scanning application will be ready to use.

### What about my patient scans on the current *eBrace* (*eBrace v4.01*)?

- All SmartSoc data will be migrated from the current *eBrace* web-portal to the new *eBrace Global* web-portal within a 48-hour period (excluding weekends) – Ref. ***eBrace Global Transition Flow***, below
- Since you can use the new configured iPhone 11 with CurveCapture O&P immediately, it is **strongly recommended to suspend any scanning activity using the old Galaxy S7(s) or iPhone 7(s)**.

#### eBrace Global Transition Flow



Notes : <sup>1</sup> 2 Hours lead-time  
<sup>2</sup> 2 Hours lead-time